

# Feverscan® Digital Reusable Thermometer

## Instructions for Use

### How to Read

1. Apply to center of forehead just above the eyebrows. Hold in place until the color stops changing, generally 15 seconds. **Read the indicator while it is still on the forehead.**
2. Green indicates the correct temperature. If green does not appear, temperature is between blue and tan.

### Intended Use

- Measure human forehead surface temperature and monitor temperature trends

### Precautions

- Environmental temperature might affect the performance of this product.

Feverscan® thermometers are to be used on the forehead area only and estimate the core body temperature by measuring the temperature of the forehead. It is an ideal way to check for changes in temperature. A Feverscan® Forehead Thermometer is not as accurate for determining the actual core body temperature as a clinical (oral / rectal) thermometer. Instead it is designed to be an unbreakable, quick reading device for temperature screening between uses of a clinical thermometer.

Forehead skin temperatures are lower than oral temperatures. Accurate to  $\pm 1^{\circ}\text{F}$  /  $\pm 0.6^{\circ}\text{C}$ , the readings on a Feverscan® Forehead Thermometer have been adjusted to approximate the difference between forehead and oral temperatures during normal to feverish conditions. Feverscan® Forehead Thermometer products indicate an approximate oral temperature. If you think that your oral temperature is different from that shown by a Feverscan® Forehead Thermometer, you should check by using a clinical thermometer. If fever is present, consult your physician for appropriate treatment.

**For External use only**

**This product contains no latex**

### Storage / Handling:

Products should be stored out of direct sunlight, dampness, or heat.

### Cleaning:

Wipe clean using soap and water or rubbing alcohol. Do not soak.

Conforms to ASTM E-1061

Manufactured by, LCR Hallcrest, Glenview, IL 60026

# Feverscan® Disposable Digital Thermometer

## Instructions for Use

### How to Read

1. Remove backing liner. Apply to center of forehead, just above the eyebrows. Hold in place until the color stops changing, generally 15 seconds. **Read the indicator while it is still on the forehead.**
2. Green indicates the correct temperature; if green does not appear, temperature will be mid way between blue and tan.

### Intended Use

- Measure human forehead surface temperature and monitor temperature trends

### Precautions

- Environmental temperatures might affect the performance of this product
- After use, discard

Feverscan® thermometers are to be used on the forehead area only and estimate the core body temperature by measuring the temperature of the forehead. It is an ideal way to check for changes in temperature. Feverscan® is not as accurate for determining the actual core body temperature as a clinical (oral / rectal) thermometer. Instead it is designed to be an unbreakable, quick reading device for temperature screening between uses of a clinical thermometer.

Forehead skin temperatures are lower than oral temperatures.

Accurate to  $\pm 1^{\circ}\text{F}$  /  $\pm 0.6^{\circ}\text{C}$ , the readings on Feverscan® thermometers have been adjusted to approximate the difference between forehead and oral temperatures during normal to feverish conditions. Feverscan® products indicate an approximate oral temperature. If you think that your oral temperature is different from that shown by Feverscan® you should check by using a clinical thermometer. If fever is present, consult your physician for appropriate treatment.

**For External use only**

**This product contains no latex**

### Storage / Handling:

Products should be stored out of direct sunlight, dampness, or heat

### Cleaning:

Wipe clean using soap and water or rubbing alcohol. Do not soak.

Conforms to ASTM E-1061

Manufactured by, LCR Hallcrest, Glenview IL 60026